

PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



Fluorescent Antifading Buffer

Cat No.: BD5014

Introduction

This product is used for protect fluorescence; it contains antifading reagent, glycerol and PBS. This product can markedly decrease the speed of light drifting process. Apply this Fluorescent Antifading Buffer as mounting reagents.

Reagents

Fluorescent Antifading Buffer 30ml

Application

For protect fluorescence

Storage & Shelf life

Store at 2-8° C for short time; -20°C for long time. Each component is stable for up to 12 Months.

NOTE

1. All the fluorescence material is easily faded with light;
2. All experiment process need keep away from light;
3. The slide protected by Fluorescent Antifading Buffer can be stored in 4°C for about 1-2 week;
4. The signal becomes decrease for long time.

Research Use

For research use only, not for use in diagnostic procedures.

MADE IN CHINA

Bioworld Technology, Inc.
1660 South Highway 100, Suite 500 St. Louis Park, MN
55416, USA. Email: info@bioworld.com
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.
No 9, weidi road Qixia District Nanjing, 210046,
P, R.China. Email: info@biogot.com
Tel: 0086-025-68037686 Fax:0086-025-68035151