



E.coli O157:H7 polyclonal antibody

Catalog: BS66115

Host: Rabbit

Reactivity: Escherichia Coli

BackGround:

Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type: H1-H56). The E. coli that cause intestinal infectious diseases including diarrhea, acute gastritis or colitis are referred to as pathogenic E. coli, which are classified into the following 4 groups according to differences in the mode of pathogenicity; enteropathogenic E. coli (EPEC), enteroinvasive E. coli (EIEC), enterotoxigenic E. coli (ETEC) and enterohemorrhagic E. coli (EHEC). Although the identification of pathogenic E. coli requires verification of their pathogenicity, pathogenic E. coli often have specific serotypes; therefore, typing of the serogroup and serotype is necessary in screening pathogenic E. coli.

Product:

0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Molecular Weight:

Swiss-Prot:

Purification&Purity:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

Applications:

ELISA:1:5000-1:10000

Storage&Stability:

Store at 4 °C short term. Aliquot and store at -20 °C long term. Avoid freeze-thaw cycles.

Specificity:

E.coli O157:H7 polyclonal antibody detects endogenous levels of E.coli O157:H7 protein.

DATA:

Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416, USA.

Email: info@bioworld.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046, P. R. China.

Email: info@biogot.com

Tel: 0086-025-68037686

Fax: 0086-025-68035151