

# **TIP47** monoclonal antibody

Catalog: MB66732

**BackGround:** 

Host:

Mouse

#### Reactivity: Human

This antibody is purified through a protein G column.

#### **Applications:**

WB (1/500 - 1/1000), FC (1/10 - 1/50)


### Storage&Stability:

Store at  $4 \,^{\circ}$  short term. Aliquot and store at  $-20 \,^{\circ}$  long term. Avoid freeze-thaw cycles.


#### **Specificity:**

Recognizes endogenous levels of TIP47 protein.

#### **DATA:**



Western blot analysis of TIP47 expression in Hela (A), K562 (B), Daudi (C) whole cell lysates.



#### Note:

For research use only, not for use in diagnostic procedure.

human placental tissue protein 17b (PP17b), binds to cytoplasmic domains of the cation-dependent (CD) and cation-independent (CI) mannose 6-phosphate receptors (MPRs) and facilitates their transport from endosomes to the Golgi complex. The inability of TIP47 to bind several proteins also transported from endosomes to the trans Golgi network indicates that TIP47 associates with a very select set of cargo molecules. In CD-MPR, TIP47 recognizes a phenylalanine/tryptophan signal sequence essential for proper sorting within the endosomal pathway.For CI-MPR binding, TIP47 requires cytoplasmic residues 48-74 of CI-MPR for high-affinity binding while residues 75-163 of CI-MPR aid in the presentation of the higher-affinity residues. Additionally, TIP47 competes with AP-2 clathrin adaptor for binding residues 24-29 of CI-MPR. In tissue extracts of cervical carcinoma patients, TIP47 is overexpressed. Dysplastic cells in high grade dysplasias express more TIP47 than dysplastic cells in low grade dysplasias, and both cytoplasmic types of dysplasias express more TIP47 than normal cervical epithelial cells. The gene encoding human TIP47 maps to chromosome 19p13.3.

Tail-interacting 47kDa protein (TIP47), known also as

### **Product:**

Mouse IgG2a kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Molecular Weight: ~ 50 kDa Swiss-Prot: O60664 Purification&Purity:

## Bioworld Technology, Inc.

 Add:
 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416,USA.

 Email:
 info@bioworlde.com

 Tel:
 6123263284

 Fax:
 6122933841

# Bioworld technology, co. Ltd.

 
 Add:
 No 9, weidi road Qixia District Nanjing, 210046, P. R. China.

 Email:
 info@biogot.com

 Tel:
 0086-025-68037686

 Fax:
 0086-025-68035151